

Middle East Technical University
Department of Physics
PHYS222 Optics and Waves Laboratory
**Experiment OW-2 Prism
Spectrometer
Laboratory Report**

Oğuzhan ÖZCAN
1852334

Partner: İnci SAİM

Teaching Assistant: Hikmet ÖZŞAHİN

March 12, 2015

Contents

1	Theory	4
2	Data and Results	6
3	Discussion and Conclusion	10
4	Application	11
5	References	12

List of Figures

1.1	Prism Spectrometer	4
1.2	Prism Dispersion	5
2.1	Refractive Index n versus Wavelength λ graph	8
2.2	Refractive Index n versus $1/\text{Wavelength}^2$ $1/\lambda^2$ graph	9
4.1	Procedure of Spectrophotometry	11

List of Tables

2.1	Data for the Prism Angle α	6
2.2	Data for the Refractive Index versus Wavelength	6
2.3	Data for wavelengths, $1/\lambda^2$ and refractive index	8

Theory

A prism spectrometer is a device where the light is dispersed by a prism. This optical element disperses parallel rays or collimated radiation into different angles from the prism. An ordinary spectrometer consists of a light source, a collimator with a slit, a prism which is located at middle of table, and an eye piece (See Figure 1.1). In this sense wavelength and index are very important because angle of deviation is depends on these quantities. The angle of deviation is shown by δ and angle of minimum daviation is shown by $\delta_m[1]$. Angle of deviation can be calculated by using

$$\delta = \theta'_2 + \omega$$

or instead of using ω we can use angle α which is formed by the edge of prism (See Figure 1.2).

$$\delta = \theta_1 + \theta'_2 + \alpha$$




Figure 1.1: Prism Spectrometer

Angle of minimum deviation δ_m is formed by the angle of incidence θ_1 . When the angle of incidence θ_1 is varied from an arbitrary point, the deviation angle δ will become smaller until it reaches a minimum point. This variation can be caused by either rotating the prism or moving the light source.




Figure 1.2: Prism Dispersion

When we make further calculations we will get some important quantities such as index of prism and relation between index n and minimum angle of deviation δ_m . For the minimum angle of deviation δ_m following equations should be satisfied;

$$n_s \sin \theta_1 = n_p \sin \theta'_1$$

$$\theta'_1 = \theta_2 = \frac{\alpha}{2} \text{ and } \theta_1 = \theta'_2 = \frac{\alpha + \delta_m}{2}$$

When the index of refraction of air to be $n_s = 1.00$

$$1.00 \sin\left[\frac{\alpha + \delta_m}{2}\right] = n_p \sin\left[\frac{\alpha}{2}\right]$$

and when we insert θ_1 and θ_2 into Snell's law, this manipulation gives a simple equation for the index of refraction of the prism [2],

$$n_p = \frac{\sin\left[\frac{\alpha + \delta_m}{2}\right]}{\sin\left[\frac{\alpha}{2}\right]}$$

In a dispersive material, the index of refraction depends on the wavelength λ of the incident light [3]. When we write the index of refraction of the prism in terms of wavelength, we get

$$n(\lambda) = \frac{\sin\left[\frac{\alpha + \delta_m(\lambda)}{2}\right]}{\sin\left[\frac{\alpha}{2}\right]}$$

Obviously the index of refraction is not same for all wavelengths or colors. As we can see from previous equation the index of refraction of prism depends on the wavelength of the incident light. We can see this relationship with *Cauchy's formula* [4],

$$n = A + \frac{B}{\lambda^2} + \frac{C}{\lambda^4}$$

where A , B and C are called Cauchy's constants.

Data and Results

Position 1	Position 2	α_i
61°28'	301°8'	59°46'
57°7'	296°21'	59°40'
65°10'	305°12'	59°39'

Table 2.1: Data for the Prism Angle α

To find the average α_i ;

$$\bar{a} = \frac{\alpha_1 + \alpha_2 + \alpha_3}{3}$$

$$\bar{a} = 59^\circ 42'$$

λ	Position 1	Position 2	δ_{min}	n
Violet	1°30'	38°52'	48°44'	1.562
Blue	1°30'	40°46'	48°35'	1.558
Green	1°30'	31°48'	48°19'	1.552
Orange	1°30'	36°43'	47°52'	1.549

Table 2.2: Data for the Refractive Index versus Wavelength

To define the related refractive index for each color we are going to use the following formula;

$$n = \frac{\sin(\frac{\alpha + \delta_{min}}{2})}{\sin \frac{\alpha}{2}}$$

for λ_{violet}

$$n = \frac{\sin(\frac{59^\circ 42' + 48^\circ 44'}{2})}{\sin \frac{59^\circ 42'}{2}}$$

$$n = \frac{\sin(\frac{108^\circ 26'}{2})}{\sin \frac{59^\circ 42'}{2}}$$

$$n = \frac{\sin(\frac{108.48^\circ}{2})}{\sin \frac{59.35^\circ}{2}}$$

$$n = \frac{\sin 54.24^\circ}{\sin 26.68^\circ} = \frac{0.811}{0.519}$$

$$n=1.562$$

We are going to have same calculation procedure for other colors.

for λ_{blue}

$$n = \frac{\sin(\frac{59^\circ 42' + 48^\circ 35'}{2})}{\sin \frac{59^\circ 42'}{2}}$$

$$n=1.558$$

for λ_{green}

$$n = \frac{\sin(\frac{59^\circ 42' + 48^\circ 19'}{2})}{\sin \frac{59^\circ 42'}{2}}$$

$$n=1.552$$

for λ_{orange}

$$n = \frac{\sin(\frac{59^\circ 42' + 47^\circ 52'}{2})}{\sin \frac{59^\circ 42'}{2}}$$

$$n=1.549$$

In this experiment we had mercury spectrum lines. Following table would be useful for rest of the calculations. I am not sure about that which wavelength should I use while plotting graph. Since we are expecting a linear slope for this graph, I am going to use values those fits to the slope.

Color	λ (nm)	$1/\lambda^2$	n
Violet	404.6	6.108×10^{-6}	1.562
Violet	407.8	6.013×10^{-6}	1.562
Blue	435.8	5.265×10^{-6}	1.558
Green	546.1	3.353×10^{-6}	1.552
Orange	577.0	3.003×10^{-6}	1.549
Orange	579.1	2.981×10^{-6}	1.549

Table 2.3: Data for wavelengths, $1/\lambda^2$ and refractive index

1. Plot a graph with indices of refraction as ordinates and wavelengths as abscissae.




Figure 2.1: Refractive Index n versus Wavelength λ graph

2. Assuming that Cauchy's equation gives the correct form of the relationship between n and λ , calculate the least square straight line for n versus $1/\lambda^2$. Plot this line on a graph with n as ordinates and $1/\lambda^2$ as abscissae. Show the experimental points on the graph. Write on the graph values A and B .




Figure 2.2: Refractive Index n versus $1/\text{Wavelength}^2 1/\lambda^2$ graph

As we can see from the graph, the graph equation is $y = 3881.3x + 1.538$. When we applied Cauchy's equation to the graph equation, A is equal to 1.538 and B is equal to 3881.3. When we solve the Cauchy's equation;

$$n = 1.538 + \frac{3881.3}{\lambda^2}$$

By solving Cauchy's equation we show that the refraction index is depend on the wavelength.

Discussion and Conclusion

1. What are the possible errors in the experiment?

The first possible error cause was the calibration of the prism spectrometer. In our setup the calibration of the prism was almost certain. The second error cause was location of the prism. Since incident ray refract from prism, location of the prism is important. In the experiment, we may had some errors. The third error cause was the reading of the Vernier Scale. It was the first time for us with using Vernier Scale. Since laboratory was dark and Vernier Scale was too small, it was very hard to read the corresponding points.

2. What kind of approximation did you take into consideration while you were obtaining the physical quantities and how do they affect your results?

Approximations those are related to experiments were the reading Vernier Scale and spectrum of light. As we know all scales have their own errors. I mean, the smallest slit was 1 mm. Besides when we look through the eye-piece, making decisions about light spectrum also caused to experimental error.

3. What discrepancies did you encounter between the calculated quantities and theoretical or literature values?

Apperantly we do not have a certain discrepancies it is because our experimental error is too small. If we calculate our percentage error by using the following formula, we will see more clearly.

$$\begin{aligned}\text{Percentage Error} &= \frac{|\text{Theoretical Value} - \text{Experimental Value}|}{\text{Theoretical Value}} \times \%100 \\ \text{Percentage Error} &= \frac{|59.42 - 59.30|}{59.42} \times \%100 \\ \text{Percentage Error} &= 0.2\%\end{aligned}$$

4. What is your overall conclusion?

We can see from the previous steps, we did experiment well. Light spectrum, wavelength and refraction index are studied well also.

Application

Spectrophotometry

Spectrophotometry is a powerful tool for analyzing substances that have color such as aspirin transition metal ion complexes and anions [5]. Principle of spectrophotometry is that every substance absorbs or transmits certain wavelengths of radiant energy but not other wavelengths. For example, chlorophyll always absorbs red and violet light, while it transmits yellow, green, and blue wavelengths. The transmitted and reflected wavelengths appear green—the color your eye “sees.” The light energy absorbed or transmitted must match exactly the energy required to cause an electronic transition (a movement of an electron from one quantum level to another) in the substance under consideration. Only certain wavelength photons satisfy this energy condition. Thus, the absorption or transmission of specific wavelengths is characteristic for a substance, and a spectral analysis serves as a “fingerprint” of the compound.




Figure 4.1: Procedure of Spectrophotometry

References

- [1] Hecht, E. (2002). *Optics* (4th ed., pp. 187-188). Reading, Mass.: Addison-Wesley.
- [2] Loshin, D. (1991). *The Geometrical Optics Workbook* (pp. 49-50). Boston: Butterworth-Heinemann.
- [3] (n.d.). Retrieved March 12, 2015, from <http://www.physics.nus.edu.sg/L1000/PC1142/PrismSpectrometer.pdf>
- [4] Chartier, G. (2005). *Introduction to Optics* (p. 363). New York: Springer.
- [5] Beran, J. (2011). *Laboratory Manual for Principles of General Chemistry* (9th ed., p. 376). Hoboken, NJ: Wiley.